

A Sensitive and Specific Test  
to Detect SARS Coronavirus

Related applications

5 The present application claims priority of U.S. Provisional Application No. 60/529,737 filed December 17, 2003.

Sequence Listing

10 The present application includes an appended Sequence Listing of 20 pages presenting 12 sequences.

Field of the Invention

15 The instant invention provides a qualitative nucleic acid amplification assay for the detection of SARS coronavirus in patient samples. The assay uses primer pairs that have been developed that provide excellent sensitivity and specificity for detection of SARS coronavirus.

20 Background of the Invention

An outbreak of atypical pneumonia, severe acute respiratory syndrome (SARS) is thought to have originated from Guang-dong Province, Republic of China in late 2002. The mortality rate of individuals suffering from SARS can be  
25 as high as 15%, depending on the age group analyzed. SARS is a highly infectious and acute condition with an extremely high mortality rate. The condition is caused by a human coronavirus, named SARS coronavirus (SARS coronavirus). The disease killed 774 patients out of 8098 probable SARS cases  
30 from November 2002 to July 2003, and has had a profound economic and social impact globally.

In many viral diseases, the spread of the virus is greatest during the early symptomatic phase that is around

BEST AVAILABLE COPY

and immediately following the onset of symptoms. Unfortunately, virus excretion is comparatively low during the initial phase of SARS. It peaks in respiratory specimens and in stools at around day 10 after the onset of the clinical illness. In order to make an early diagnosis, it is therefore necessary to use highly sensitive tests that are able to detect the low levels of viral genome present during the first days of the illness.

There are many non-standardized and sensitive tests under development in many countries. The available SARS RT-PCR based diagnostic tests often suffer the drawback of being complex and difficult to administer. The typical SARS diagnostic test uses nested (two step) polymerase chain reaction (PCR) to accomplish a certain level of specificity and sensitivity. See, e.g., "SARS-CoV Specific RT-PCR Primers", by William J. Bellini, Ph.D. Chief, Measles Virus, Section DVRD/NC1D/CDC, CDCprimers.pdf, obtainable from the World Health Organization (WHO), which is hereby incorporated by reference in its entirety, for description of the typical PCR test for SARS.

#### Brief Description of the Drawings

Figure 1 shows the portions of the SARS coronavirus genome amplified by the IMCB primer sets.

Figure 2 shows the portion of the SARS coronavirus genome amplified by the IMCB-3 primer set and aligns the IMCB-3 primers and the IMCB-3 probe along the sequence of the SARS coronavirus genome. The upper strand sequence is shown as nucleotides 4609-4765 of SEQ ID NO: 1. The lower strand is shown as SEQ ID NO: 12.

Figures 3A-3C show gels demonstrating the efficacy of the primers of the invention. Figure 3A shows the detection of SARS coronavirus using the IMCB-2 primer set where the

virus copy number per sample loaded varies from 26.1 copies to 0.07 copies. Lane 1 is a marker, lanes 2 & 3 contain 26.1 copies of the virus per 5  $\mu$ l, lanes 4 & 5 contain 12.6 copies of the virus per 5  $\mu$ l, lanes 6 & 7 contain 1.96 copies of the virus per 5  $\mu$ l, lanes 8 & 9 contain 2.0 copies of the virus per 5  $\mu$ l, lanes 10 & 11 contain 0.07 copies of the virus per 5  $\mu$ l, lane 12 contains a negative control of an unrelated virus and lane 13 contains another marker.

Figure 3B shows a second experiment detecting SARS coronavirus using the IMCB-2 primer set, where the virus copy number per sample loaded varies from 26.1 copies to 0.08 copies. Lane 1 is a marker, lanes 2 & 3 contain 26.1 copies of the virus per 5  $\mu$ l, lanes 4 & 5 contain 8.2 copies of the virus per 5  $\mu$ l, lanes 6 & 7 contain 2.6 copies of the virus per 5  $\mu$ l, lane 8 contains 0.8 copies of the virus per 5  $\mu$ l, lane 9 contains 0.25 copies of the virus per 5  $\mu$ l, lanes 10 & 11 contain 0.08 copies of the virus per 5  $\mu$ l, lane 12 contains a negative control of an unrelated virus and lane 13 contains another marker.

Figure 3C shows the amplified product from a sample containing 5 copies of SARS coronavirus genomic RNA per run in a total of 5  $\mu$ l (duplicated). The product is resolved by 3% agarose gel electrophoresis. A 10% of total reaction volume (5  $\mu$ l) is loaded per lane. Lane 1, amplified product; lane 2, amplified product of a duplicate reaction; lane M, 100 bp ladder.

Figures 4A to 4C show the sensitivity achieved using the present invention to detect SARS coronavirus nucleic acid with the primer set IMCB-1. Figure 4A shows results achieved with 8.8 pfu (2200 copies) per sample (lanes 1 and 2) to 0.08 pfu (22 copies) per sample (lanes 5-6). Lanes 7 and 8 show a no virus control. Figure 4B shows the results of another run of the same assay using from 0.08 pfu (22

copies) per sample to 0.0008 pfu (0.2 copies) per sample. Lanes 11 and 12 show a no virus control. Figure 4C shows a third run using from 0.08 pfu (22 copies) per sample to 0.004 pfu (1 copy) per sample. Lanes 7 and 8 are a no virus control sample. M is a molecular length marker.

### Description of the Invention

At the time of the original SARS outbreak there was a lack of rapid detection. Sensitive and specific rapid detection would have allowed quick diagnosis of infected patients to enable better containment of the spread of the epidemic. A PCR-based assay was developed at the Bernhard Nocht Institute (BNI) (Drosten et al, 2003) and is available commercially from Artus. The primers identified by BNI and used by Artus are from the SARS coronavirus non-structural protein 9 that encodes an RNA polymerase. According to WHO recommendations, results of these tests should still not be used to rule out a suspected case of SARS (WHO Update 71).

Because presently available tests are not generally able to detect the requisite small amounts of SARS coronavirus (SARS coronavirus), they do not yet play a role in patient management and case control, as SARS patients may be capable of infecting others during the initial phase and therefore need to be reliably detected and quickly isolated.

Coronaviruses are a family of RNA viruses with a large envelope that propagate in the cytoplasm of host cells and usually cause mild respiratory disease in man and animals.

The SARS Coronavirus has been isolated and sequenced. A prototype sequence of 29,727 basepairs can be found at GENBANK, under Accession No. AY278741, hereby incorporated by reference and presented also as SEQ ID NO: 1. See also, Y.J. Ruan et al., *Lancet* 361:1779-1785 (2003), analyzing the genome sequence of 14 different isolates, and P.A. and Rota

et al., *Science* 300:1377-1378 (2003), characterizing one of the first isolates to be associated with SARS.

Sequencing of the complete genome of the SARS virus from a number of different isolates has indicated that the virus has a typical coronavirus genome organization, but that the virus is not closely related to any other known coronaviruses.

The SARS virus encodes 14 open reading frames (ORFs), including the replicase 1a and 1b proteins and four structural proteins, spike protein (S), envelope protein (E), membrane protein (M) and nucleocapsid protein (N).

Ruan et al. (2003, *supra*) compared the genome sequence of a Singapore case of SARS coronavirus and a database of other coronavirus genomic sequences. From this they were able to find which regions of SARS coronavirus were homologous to other coronavirus and thus conserved among coronavirus strains and which sequences were unique to SARS coronavirus Singapore strain SIN2500.

The instant invention provides a simple, sensitive and specific diagnostic test. This test provides a simple yet sensitive and specific nucleic acid amplification system compared with others that have so far been developed. By use of the primers described herein, the instant invention is made more sensitive and specific than the detection methods of the prior art. Furthermore, such specificity and sensitivity may be enhanced by using a one step PCR method, instead of a two step PCR.

The present invention utilizes specific primer pairs designed from the SARS coronavirus non-structural protein 1 (NSP1) a putative proteinase. These primers can be used with many techniques to detect the presence of SARS coronavirus.

In one embodiment the instant invention provides a simple gel-based RT-PCR detection kit. Such a kit will

include one or more primers and/or probes according to the invention, for example a kit may contain primers consisting of one or more polynucleotides comprising a nucleotide sequence of SEQ ID NOS: 3, 4, 6, 7, 9, 10 and 11. A kit  
5 according to the invention may optionally include a positive control nucleic acid, for example a SARS coronavirus genomic nucleic acid, or at least a portion thereof comprising the NSP1 region, as either RNA or DNA.

The present invention also provides a method for detecting  
10 SARS coronavirus nucleic acid in a sample. The method may be generally described as comprising amplifying a nucleic acid of a sample with a reverse transcriptase and at least one primer specific for the NSP1 region of a SARS coronavirus to generate a nucleic acid amplification product.

15 The amplification product is then analyzed, seeking to detect an expected nucleic acid amplification product. Detection of the expected product indicates the presence of SARS coronavirus nucleic acid in the sample.

The primers of the instant invention can also be used as a  
20 primer set for real-time PCR detection using PCR platforms such as the Roche LightCycler™, the Stratagene Real-time PCR system, the Applied Biosystems ABI 7000 real time PCR analyzer or any other suitable detection platform.

Specific primer pairs have been designed by comparing  
25 conserved regions among SARS coronavirus strains and avoiding sequence regions that were conserved among all coronaviruses generally. A portion of the SARS coronavirus genome that encodes the NSP1 proteinase was selected by this comparison as a detection target.

As illustrated in Figure 1, the primers of the present invention hybridize in the portion of the SARS coronavirus genome from about nucleotide number 2200 to about nucleotide 9800 in a manner as to amplify this region, or a portion

thereof. This region is known as the NSP1 region, which codes for a putative proteinase. This region was chosen as a putative region to amplify because unlike the RNA polymerase region utilized in prior art assays, there is a significant portion of it in which no mutations appear to have occurred among SARS coronaviruses and this region is very specific to all isolates of SARS coronavirus; a hypothesis which has been examined by an NCBI (National Center for Biotechnology Information) Blast search. In the region from about nucleotide 2650 to about nucleotide 7850 has been identified that appears to bear no mutations among several strains (Ruan et al., 2003). A preferred part of the NSP1 region for amplification is the part from and including nucleotide 4609 to and including nucleotide 7003.

The primers of the invention should be at least 16 nucleotides in length, more preferably at least 18 nucleotides in length, still more preferably at least 20 nucleotides in length. The primers should be less than 50, preferably less than 30, more preferably less than 25 nucleotides in length, so as to preserve the specificity of the primers.

Thus, the SARS detection method of the invention lies generally in use of a set of primers that are specific for the NSP1 region of the SARS coronavirus genome for PCR amplification of this part of the genome and detection of the amplification product. The method of the invention can be performed, for example by amplifying nucleic acids present in the sample using a forward primer and a reverse primer selective for the region of the SARS genome from nucleotide 6652 to nucleotide 7003, or using a forward primer and a reverse primer selective for the region of the SARS genome from nucleotide 4609 to nucleotide 4765, said primers having a certain primer length in nucleotides and

being separated by a separation length that is a certain number of nucleotides, to obtain an amplification product. The amplification product is then detected. The amplification product can be detected, for example, by  
5 determining the length of the amplification product in nucleotides, either by a chromatographic method or by a gel electrophoretic method, e.g by electrophoresis in 2 or 3% agarose. The presence of an amplification product having a  
10 length in nucleotides that is the sum of the forward primer length, the reverse primer length and the separation length indicates the presence of SARS coronavirus nucleic acid in the sample.

Alternatively, the product can be detected using a hybridization probe, for example using real-time fluorescent detection in the Taqman™ system. The hybridization probe preferably comprises a nucleotide sequence that is the same as that of a portion of the amplification product that would be obtained using the amplification primers selected and a SARS coronavirus genomic nucleic acid as a template. The hybridization probe should be at least 16 nucleotides in length, more preferably at least 18 nucleotides in length, still more preferably at least 20 nucleotides in length. The probe should be less than 50, preferably less than 30, more preferably less than 25 nucleotides in length, so as to preserve the specificity of the probe.

The essential function of a primer or probe according to the present invention is to specifically hybridize to a SARS coronavirus nucleic acid, either an RNA or DNA, and not to cross-hybridize to other coronavirus nucleic acids or to nucleic acids of other viruses. Thus, a primer or probe according to the present invention "consists essentially of" a nucleotide sequence if it includes that sequence and additional nucleotides that do not impair the ability of the



primer or probe to specifically hybridize to a SARS coronavirus nucleic acid under the conditions selected for performing a diagnostic assay according to the invention.

#### Materials and methods generally used in the Examples

##### **SARS coronavirus culture**

The SARS coronavirus isolate (2003VA2774), which has been previously sequenced (Ruan et al, 2003), is used for this study. The virus stock is propagated in Vero E6 cells (ATCC: C1008) with medium 199 (Sigma Aldrich, USA) supplemented with 5% fetal calf serum (FCS) (Biological Industries, Israel). When more than 75% of the cell monolayer showed cytopathic effects, the culture supernatant is harvested, clarified by centrifugation at 1300 x g, aliquoted and stored at -80°C until use. The PFU of the current preparation is determined as  $1 \times 10^7$  PFU/ml.

##### **Plaque assay**

A plaque assay is carried out to determine the virus titer in the culture. Briefly, 100 µl series of 10-fold dilution of the virus stock is added to a confluent monolayer of Vero E6 cells in a 24-well plate and incubated for 1 hour at 37°C.

Following this, 1 ml of 1% carboxymethylcellulose overlay in medium 199 with 5% FCS is then added to each well. After 4 days of incubation at 37°C in 5%CO<sub>2</sub>, the cells are then fixed with 10% formalin and stained with 2% crystal violet. The plaques are counted visually and the virus titer determined.

##### **RNA extraction**

##### SARS coronavirus Standards:

10-fold dilution of the stock virus is prepared in serum obtained from a healthy volunteer. RNA is extracted

using the QIAGEN Viral RNA Kit (QIAGEN GmbH, Germany) according to the instructions given in the product insert.

Patient specimens:

Virus isolation is performed on a bronchoalveolar lavage specimen of SARS cases belonging to the original case cluster from Singapore. RNA is directly extracted from the specimen using a Qiagen QIAamp viral RNA extraction kit (catalog no. 52906) according to the instructions given in the product insert.

10 Other Viruses:

RNA is directly extracted from the stock vial obtained from ATCC (VA, USA) using the QIAGEN Viral RNA Mini Kit (QIAGEN GmbH, Germany) according to the instructions given in the product insert.

15 MRC-5 cell line:

Total RNA is extracted directly from the normal diploid human fibroblast cell line MRC-5 (ATCC CCL171) using a Qiagen RNA extraction kit (catalog no. 74104) and RNA is quantitated using a spectrophotometer.

20

Example 1: Primer design

Primers are selected using the sequence of the SARS coronavirus Urbani strain (AY278741) based on conserved regions of SARS coronavirus genome (Rota et al. 2003).

A set of primer pairs in the proteinase gene region (position 6652-7003, of a non-structural protein 1 (NSP1) region) is found to be most suitable as it exhibits the lowest cross homology with other viruses (Ruan et al. 2003). The primers are designed to take into account possible mismatches throughout the genome and to avoid or at least minimize primer dimer formation.

The NSP1 region (proteinase) target of the primers is generally well conserved among isolates of SARS coronavirus.

In most of the NSP1 region, there is no match between a SARS fragment and other coronaviruses compared with NSP9 (RNA polymerase), which has very strong homology across the coronavirus group (Ruan et al. 2003).

5 In addition, no mutation has been found occurring in a significant portion of this region (Ruan et al. 2003). It is also very specific to SARS coronavirus (all known isolates), as confirmed by testing amplification of RNA samples from other viruses using primers targeting this  
10 region (Table 1) and by NSTBI BLAST search of all known sequences in the NSTBI database excluding SARS coronavirus sequences.

Three sets of primer pairs are identified from the NSP1 proteinase region. The location of these three sets in  
15 relation to the SARS coronavirus genome can be seen in Figure 1. The first set amplifies a segment IMCB-1 that is 352 base pairs in length, SEQ ID NO: 2. This sequence is flanked by an upper primer (IMCB-1-U, SEQ ID NO: 3) and a lower primer (IMCB-1-L, SEQ ID NO: 4). These primers can be  
20 used to specifically detect the presence of SARS coronavirus nucleic acids in a sample.

IMCB-1-U (19-mer) : 5'ACATCAAATTGCGCTAAGA3'

(SEQ ID NO: 3)

25

IMCB-1-L (21-mer) : 5'ACAATTCTCTAACGCCATTAC3'

(SEQ ID NO: 4)

Set 2 relates to a fragment called IMCB-2 of 157 base  
30 pairs in length, SEQ ID NO: 5. This sequence is flanked by an upper primer (IMCB-2-U, SEQ ID NO: 6) and a lower primer (IMCB-2-L, SEQ ID NO: 7). These primers can be used in like manner to the IMCB-1 primer set to specifically detect

the presence of SARS coronavirus nucleic acids in a sample. Both the IMCB-1 and IMCB-2 primer sets can be used with the reagents and conditions set forth, for instance, in the One-Step RT-PCR described in Example 3.

5

IMCB-2-U (19-mer): 5'GCCGTAGTGTCTAGTATCAT3'

(SEQ ID NO: 6)

IMCB-2-L (21-mer): 5'CACCTAACTCTGTACGCTGTC3'

(SEQ ID NO: 7)

10

Set 3 relates to a fragment called IMCB-3 that is a portion of IMCB-2 that is 77 base pairs in length, SEQ ID NO: 8. This sequence is flanked by an upper primer (IMCB-3-U, SEQ ID NO: 9) and a lower primer (IMCB-3-L, SEQ ID NO: 10). A probe oligonucleotide (IMCB-3-probe, SEQ ID NO: 11) is situated between the two primers of IMCB-3. The primers and the probe are used to detect the presence of SARS coronavirus, for example in a real-time assay using fluorescent detection of the amplified fragment. The lower primer of set 2, IMCB-2-L and of set 3, IMCB-3-L, are almost identical, but IMCB-3-L is 3 nucleotides longer at the 3' end. These primers and/or probe can be used to specifically detect the presence of SARS coronavirus nucleic acids in a sample.

25

IMCB-3-U (26-mer): 5'GCACTTTGTAGAAACAGTTTCTTTGG3'

(SEQ ID NO: 9)

IMCB-3-L (24-mer): 5'CACCTAACTCTGTACGCTGTCCTG3'

(SEQ ID NO: 10)

30

IMCB-3-Probe (17-mer): 5'TGGCTCTTACAGAGATT3'

(SEQ ID NO: 11)

Example 2: Detection of SARS coronavirus using the IMCB-1 primer pair

The first primer pair (IMCB-1) described above is designed for amplifying the portion of the SARS coronavirus genome from nucleotide 6652 to nucleotide 7003. The forward primer, also herein called an "upper primer" is a 19-mer that hybridizes to the SARS coronavirus genome beginning at position 6652 and has the sequence of SEQ ID NO: 3.

The reverse primer, called herein a "lower primer" is a 21-mer that hybridizes to the SARS coronavirus genome at position 6983 and has the sequence of SEQ ID NO: 4.

RNA is extracted from samples thought to contain SARS RNA by known methods. The RNA is then converted to DNA using a reverse transcriptase or any other method known in the art. A sample mixture is converted into cDNA in a typical manner using a 1st Strand cDNA Synthesis Kit for RT-PCR (sold by Roche, Basel, Switzerland, catalog no. 1 483 188).

The first strand cDNA reaction is carried out using the following reagents at the indicated concentrations:

<u>REAGENTS</u>	<u>FINAL CONC.</u>
1. 10x Reaction buffer	1x
2. MgCl <sub>2</sub> , 25 mM	5 mM
3. dNTP mix, 10 mM ea.	1 mM ea.
4. Specific primer, 20 µM	1 µM (IMCB-1 Lower Primer)
5. RNase inhibitor	50 units
6. AMV reverse transcriptase	20 units
7. Gelatin (0.5 mg/ml)	0.01 mg/ml
8. Sterile water	-
9. RNA sample	-

to produce a total reaction volume of 20.0 µl. The reaction is allowed to proceed at 25°C for 10 minutes and then at 42°C for 60 minutes.

The reverse transcriptase is then inactivated by incubating the reaction at 99°C for 5 minutes and cooling to 4°C for 5 minutes.

The cDNA is then amplified by adding 1 µl of 1 µg/µl of single stranded DNA to each sample and preparing the sample for PCR. It will be recognized by those of skill in the art that other methods of amplification known in the art can be performed to amplify the DNA.

PCR is performed using the following reagents and conditions. The reaction mixture is then prepared using the following reagents and concentrations.

<u>REAGENTS</u>	<u>FINAL CONC.</u>
1. 2x master mix*	x1
2. Upper primer	0.3 µM
3. Lower primer	0.3 µM
9. DNA sample	-

Total reaction volume            50.0 µl

\*Promega #M7501 (Madison, Wisconsin)

The thermal cycling is performed using a Stratagene (La Jolla, California) Robocycler 96 for the respectively enumerated steps, temperatures and times.

1. Initial denature,    95°C 5 min
2. Denaturation,        95°C 45 sec
3. Annealing,            49°C 80 sec
4. Extension,            72°C 50 sec

= 2,3,4 cycling: 35=

5. Final extension,    72°C 3 min

The PCR product is analyzed on a 2.0% agarose gel. The  
5 detection can be alternatively done by an ABI (Foster City, California) PRISM 7000 Sequence Detection System to confirm the presence of the correct amplified region, which is identified as a nucleic acid fragment of 352 nucleotides.

It will be obvious to a skill person in the art that the IMCB-1 primer sets can be used in the Two-Step RT-PCT as well as in the One-Step RT-PCR. In particular, the IMCB-1 primer sets may be used with the reagents and conditions set forth, for instance, in the One-Step RT-PCR described in Example 3.

Example 3: SARS coronavirus detection using IMCB-2 primer set

The second primer pair set described above is designed for amplifying the portion of the SARS coronavirus genome from nucleotide 4609 to nucleotide 4765. The forward primer, called herein "Upper Primer" or "IMCB-2-U" is a 19-mer that hybridizes to the SARS coronavirus genome beginning at position 4609 has the following sequence:

IMCB-2-U 5' GCCGTAGTGTCAGTATCAT3' (SEQ ID NO: 6)

The reverse primer, called herein "Primer Lower" or "IMCB-2-L" is a 21-mer that hybridizes to the SARS coronavirus genome at position 4765 has the following sequence:

IMCB-2-L 5' CACCTAACTCTGTACGCTGTC3' (SEQ ID NO: 7)

The method to perform the assay is as follows.

5  $\mu$ l of RNA sample is diluted to a 50  $\mu$ l reaction volume with 45  $\mu$ l of a premixture solution containing reaction buffer, Q-solution (Qiagen, catalog no. 210210) and dNTP mix, at a final concentration of 400  $\mu$ M each, and an upper primer and a lower primer at a final concentration of 0.6  $\mu$ M each. RNase inhibitor at 10 units/ reaction, enzyme mix, RNase free water is also added.

The thermal cycling is performed using a Stratagene Robocycler 40 (La Jolla, California) with the following

steps, reverse transcription at 50°C for 30 min and initial denature at 95° C for 15 min, followed by denaturation at 95° C for 45 sec, annealing at 50°C for 80 sec, extension at 72°C for 50 sec. The cycle is repeated 42 times. The entire RT PCR is completed after a final extension at 72° C for 10 min. The PCR product is analyzed by gel electrophoresis in 2.0% agarose.

With patients' samples and pre-determined SARS coronavirus standards, it is found that the detection limit of the diagnostic test of the instant invention is approximately 200 copies/ml (1 copy/ 5 µl reaction) for the virus as confirmed when measured using ARTUS RealArt™ HPA-Coronavirus LC RT PCR Reagents (cat No: 5601-03).

#### 15 Example 4: A RT-PCR SARS diagnostic kit and its use

A kit according to this Example is typically prepared to contain 50 or 100 reactions. The kit is composed of the items listed below; it should be stored at -20°C in a non-frost-free freezer. Figure 3C shows results that are obtained using this kit with the IMCB-2 primer set.

The Reverse Transcriptase-Polymerase Chain Reaction (RT-PCR) kit described herein is used for detecting the presence of severe acute respiratory syndrome coronavirus (SARS coronavirus) RNA in samples extracted from specimens with an appropriate RNA extraction method of choice. The region of the SARS coronavirus genome amplified in this assay lies in the the NSP-1 region (proteinase) of the virus genome. This kit is optimized to detect a few molecules of the viral RNA in 5 ul of test sample and the entire procedure is performed in one step.



*Components*

- The kit of this example consists of the following 4 tubes:

<u>Tube No.</u>	<u>Component</u>	<u>Storage Conditions</u>
(Short, Long term)		
5 Tube 1	RT-PCR Enzyme Mix (e.g. from Qiagen)	-20°C, -20°C
10 Tube 2	RT-PCR Reaction Mix (e.g. from Qiagen)	4°C , -20°C
Tube 3	Primer Mix (Upper & Lower) (30 µM each primer in 20 mM Tris, 1 mM EDTA pH 8.2)	4°C , -20°C
15 Tube 4	Positive Control (RNA transcripts of the gene targeted by the primers)	-20°C, -80°C

*Protocol, One Step RT-PCR**1. Sample Preparation*

- 20 In a RNase-Free Eppendorf tube (0.5 ml or 0.2 ml size),  
add the following reagents per test/per reaction:

	<u>Description</u>	<u>50 µl/Rxn</u>	<u>20 µl/Rxn</u>
25	RT-PCR Enzyme Mix	2.0 µl	0.8 µl
	RT-PCR Reaction Mix	42.0 µl	16.8 µl
	Primer Mix (U & L)	1.0 µl	0.4 µl
	RNA Samples	5.0 µl	2.0 µl
	Total Volume	50.0 µl	20.0 µl

*30 2. Thermal cycling protocol-A*

This thermal cycling protocol is used for three-block  
type PCR cyclers such as the RoboCycler® by Stratagene:

<u>Stage</u>	<u>Temp (°C)</u>	<u>Duration</u>	<u>No.</u> <u>Of</u> <u>Cycle(s)</u>	<u>Step</u>
1	50	30 mins	1	Reverse transcription
5 2	95	15 mins	1	Initial denaturation
3	95	45 secs	42	Denaturation
	57	80 secs		Annealing
	72	50 secs		Extention
4	72	10 mins	1	Final Extention

10

### 3. Thermal cycling protocol-B

This thermal cycling condition is for one-block type PCR cycler such as the Px2 Thermal Cycler by Thermo Electron.

<u>Stage</u>	<u>Temp (°C)</u>	<u>Duration</u>	<u>No.</u> <u>Of</u> <u>Cycle(s)</u>	<u>Step</u>
1	50	30 mins	1	Reverse transcription
2	95	15 mins	1	Initial denaturation
20 3	95	18 secs	42	Denaturation
	57	36 secs		Annealing
	72	33 secs		Extention
4	72	10 mins	1	Final Extention

25 It is preferred to put mineral oil in the wells of the thermocycler to maximize the conduction of the heat between a reaction tube and a well if it is necessary.

### 4. Termination of PCR reaction

30 This step is optional.

- (1) Add 30 µl of chloroform/tube. Vortex mix for 5 seconds.

- (2) Centrifuge for 2 minutes. (Top=Aqueous phase,  
Bottom=Organic phase)

The sample is reserved as the top, aqueous phase.

5 *5. Detection by electrophoresis*

The products of the PCR reactions are resolved by DNA gel electrophoresis by using 5 µl of the product reaction mixture product per lane. 3% agarose gel provides good resolution; gels are typically run at 100 V for 30 min.

10 The expected product size is 157bp when the IMCB-2 primer set of Example 1 is used. A result of such an assay is shown in figure 3C. This example shows the amplified product from the sample containing 5 copies /rxn (5 µl) SARS coronavirus RNA run in duplicate. The product was resolved  
15 by 3% agarose gel electrophoresis. 10% of the total reaction volume (5 µl) was loaded per lane. Lane 1, product of run 1; Lane 2, product of duplicate run; lane M, 100 bp ladder.

20 Example 5: Specificity of RT-PCR using the IMCB primers

To verify that the primer sets designed in Example 1 can be used to detect SARS coronavirus specifically; the amplification of selected viruses is tested by RT-PCR using the IMCB primer sets 1 and 2. The following viruses are  
25 tested at the indicated titer to check the specificity of the IMCB RT-PCR primer sets 1 and 2:

Human Coronavirus 229E, ATCC VR-740	2.8 x 10 <sup>(6)</sup> PFU/ml
Human Coronavirus OC43, ATCC VR-759	3.5 x 10 <sup>(7)</sup> LD <sub>50</sub> /0.02 ml

Avian infectious bronchitis virus	5.8 x 10(6) PFU/ml
NCBI M95169	
Dengue virus, NCBI M87512	1.1 x 10(6) PFU/ml
5 Yellow Fever virus, vaccine strain 17D	0.5 x 10(6) PFU/ml
Human Enteric Coronavirus, ATCC VR-1475	0.2 x 10(6) PFU/ml
Bovine Coronavirus, ATCC VR-874	1.1 x 10(6) PFU/ml
Rabbit Coronavirus, ATCC VR-920	23.3 x 10(6) PFU/ml
Mouse hepatitis virus, ATCC VR-764	0.6 x 10(6) PFU/ml
10 Canine Coronavirus, ATCC VR-809	36.8 x 10(6) PFU/ml
Rat Coronavirus, ATCC VR-1410	1.1 x 10(6) PFU/ml
Feline-CoV RNA, ATCC VR-989	1 x 10(5.5) CID[50]/ml

The results in Table 1 show the IMCB primer pair specificity. IMCB primer pair 1 and IMCB primer pair 2 together with the ARTUS detection methods are used for detection of other Coronaviruses and viruses unrelated to SARS coronavirus. Negative results are observed for all the viruses using the IMCB-1 and IMCB-2 primer pairs. Thus, it is demonstrated that the IMCB pairs are highly specific to SARS coronavirus for detection.

Table 1: Specificity of IMCB-1 and IMCB-2 primer pairs

Sample	Conc.	1st ARTUS	2nd ARTUS	IMCB 1	IMCB 2
MRC5 RNA Ex1	2 ug/ul	INCONCLUSIVE	Neg	Neg	Neg
MRC5 RNA Ex2	1 ug/ul	Neg	na	Neg	Neg
Avian infectious bronchitis virus, NCBI M95169	5.8 x 10(6) PFU/ml	Neg	na	Neg	Neg
Dengue virus, NCBI M87512	1.1 x 10(6) PFU/ml	Neg	na	Neg	Neg
Yellow Fever virus, vaccine strain 17D	0.5 x 10(6) PFU/ml	Neg	na	Neg	Neg
Human Enteric coronavirus, ATCC VR-1475	0.2 x 10(6) PFU/ml	Neg	na	Neg	Neg
Bovine coronavirus, ATCC VR-874	1.1 x 10(6) PFU/ml	Neg	na	Neg	Neg
Rabbit coronavirus, ATCC VR-920	23.3 x 10(6) PFU/ml	Neg	na	Neg	Neg
Mouse hepatitis virus, ATCC VR-764	0.6 x 10(6) PFU/ml	Neg	na	Neg	Neg
Canine coronavirus, ATCC VR-809	36.8 x 10(6) PFU/ml	Neg	na	Neg	Neg
Rat coronavirus, ATCC VR-1410	1.1 x 10(6) PFU/ml	Neg	na	Neg	Neg
Human coronavirus 229E, ATCC VR-740	2.8 x 10(6) PFU/ml	Neg	na	Neg	Neg

Example 6: Comparison of the sensitivity of the available  
SARS coronavirus detection kits

The available SARS detection kits available, Eiken, Artus and Roche SARS diagnostics, are compared with the three IMCB primer pairs in their ability to detect the same virus standard. The Eiken kit is tested using a one step RT-Lamp. (T. Notomi et al., "Loop-mediated Isothermal Amplification of DNA", *Nucleic Acids Research* 15:E63 (2000).) The IMCB-1 primer pairs are tested using the concentrations of reagents set forth in Example 2. However the reagent materials are those described in Example 3 and the PCR protocol was that of the one step RT-PCR described in Example 3. The IMCB-2 primer pairs are tested using the reagent materials and concentrations and the one step RT-PCR protocol described in Example 3. The amplification products that are obtained using the IMCB-1 and IMCB-2 primer sets are analysed by agarose gel electrophoresis with ethidium bromide staining. The IMCB-3 primer pair is tested using a one step RT-PCR kit prototype optimized for the ABI 7000 Real Time system using a Taqman<sup>TM</sup> probe described in Example 1. The Artus ABI kit is tested using RealArt<sup>TM</sup> HPA-Coronavirus TM RT PCR, Abbot List No. B3K360 REV.2003-10-R2. (10 µl sample is used per reaction instead of 5 µl). The Artus light cycler kit is tested using RealArt<sup>TM</sup> HPA-Coronavirus LC RT PCR Reagent (Cat No: 5601-03). The Roche kit is tested using Light Cycler<sup>TM</sup> SARS Quantification kit (cat 03604438001) Version 1.

Table 2 shows the results of the comparative testing. The viral sample copy numbers used is decreased from 83 copies per 5 µl at the top of the table to 0.1 copies per 5 µl at the bottom of the table in column 1 under the sample volume/run (vol/rxn). Across the table the numbers in

brackets behind Pos indicate how many positive detections there are per the total number of tests done.

From Table 2 it is seen that the IMCB-2 primer sets provide the most sensitive detection; the IMCB-1 and -2 primer pairs can reliably detect as little as 0.8 copies of viral nucleic acid per 5  $\mu$ l sample.

Table 2

Detection of SARS-CoV by BKEN, IMCB, ARTUS, and ROCHE Systems												
	BKEN LAMP		IMCB 1	IMCB 2	IMCB 3 (ABI)			ARTUS (ABI)				ROCHE
Sample code Sample vol/rxn	Run1 5 ul	Run2 5 ul	Run 1 5 ul	Run 1 5 ul	Run1 5 ul	Run2 5 ul	Run 1 10 ul	Run 1 5 ul	Run 2 5 ul	Average (copy/ul)	STDV, CV	Run 1 5 ul
(N-5.5) 83 copies/5ul	Pos(8/8)	Pos(6/6)	Pos (12/12)	Pos (8/8)				13.080	13.460	12.164	0.923	
								12.970	11.740	Pos (7/7)	7.6%	
								10.840	11.360			
								11.700				
(N-6.0) 26 copies/5 ul	Pos(8/8)	Pos(6/6)	Pos (12/12)	Pos (8/8)	Pos (12/12)	Pos (6/6)		5.236	6.114	5.215	0.457	Pos(8/8)
								5.327	4.534	Pos (7/7)	8.8%	
								5.192	4.806			
								5.296				
(N-6.5) 8 copies/5 ul	(not done)	Pos(4/6)	Pos (12/12)	Pos (8/8)	Pos (12/12)	Pos (6/6)		2.894	2.094	2.520	0.479	(not done)
								3.335	2.780	Pos (7/7)	19.0%	
								2.491	2.202			
								1.847				
(N-7.0) 3 copies/5 ul	Pos(4/8)	Pos(4/6)	Pos (12/12)	Pos (8/8)	Pos (12/12)	Pos (6/6)	Pos (2/2)	0.354	0.472	0.393	0.232	Pos(0/8)
								0.672	0.192	Pos (7/7)	59.0%	
								0.700	0.549			
								0.013	0.191			
(N-7.5) 0.8 copies/5 ul		Pos(0/6)	Pos (12/12)	Pos (8/8)	Pos (10/12)	Pos (6/6)	Pos (2/2)	0.417	0.219	0.413	0.318	
								0.000	0.982	Pos(6/7)	77.0%	
								0.758	0.192			
								0.324				
(N-8.0) 0.3 copies/5ul	Pos(0/8)		Pos (4/10)	Pos (1/8)	Pos (2/12)	Pos (1/6)	Pos (0/2)	0.000	0.000	0.014	0.034	
							noise*	0.000	0.097	Pos(1/7)	244.9%	
								0.000	0.000			
								0.000				
(N-8.5) 0.1 copies/5 ul			w Pos (2/10)	F(2/6)	Pos (1/12)	Pos (1/6)		0.000	0.000	0.061	0.137	
								0.368	0.000	Pos(1/6)	223.6%	
								0.000	0.000			
*background noise moderately high												
BKEN: One step RT-LAMP												
IMCB 1: One step RT-PCR with Primer set #1, product was analysed by Agarose gel/EBr.												
IMCB 2: One step RT-PCR with Primer set #2, product was analysed by Agarose gel/EBr.												
IMCB 3 (ABI): One step RT-PCR kit prototype optimized for ABI 7000 Real Time System using TaqMan Probe												
ARTUS (ABI) : RealArt HPA-Coronavirus TM RT PCR, ABBOT Ltd.No. B3K360 REV.2003-10-R2. Note: 10 ul sample was used per reaction instead of 5 ul.												
ARTUS: RealArt HPA-Coronavirus LC RT PCR Reagents (cat No: 5601-03)												
ROCHE Light Cycler SARS Quantification Kit (cat 03604438001) Version 1												

Example 7: Analysis of patient samples

Clinical samples are obtained from a number of patients and are analyzed by the assay of the invention, using primer set IMCB-2 and the assay method described in Example 3. The  
5 assay method described in Example 4 is used for samples analysed by Artus using their PCR kit. Results are shown in Table 3.

Column 1 in Table 3 is a sample identification number. Column 2 is the description of the type of sample that was  
10 taken from the patients. Columns 3 and 4 indicate the results of the Artus detection method. Column 5 is the results of the detection using RT PCR with IMCB-2 primer set. The final column 6 documents notes referring to each sample.

Table 3: Analysis of patient samples

INCB ID	DESC	1ST ARTUS	2ND ARTUS	PS 2 INCB	INCB Note
P65	PLASMA	-N		N	
P66	STOOL	N		N	
P67	STOOL	247		N	
P68	STOOL	-N		N	
P69	STOOL	194		N	
P70	STOOL	N		N	
P71	STOOL	-N		POS	a->g.(4656)
P72	STOOL	N		N	
P73	STOOL	N		N	
P74	STOOL	N		N	
P75	STOOL	N		N	
P76	THROAT SWAB	N		N	
P77	STOOL	N		N	
P78	BAL	-N		N	
P79	SALIVA	N		N	
P80	SERUM	N		N	
P81	PLASMA	N		N	
P82	STOOL	INC		N	
P83	SPUTUM	N		N	
P84	SPUTUM	N		N	
P85	STOOL	N		N	
P86	SALIVA	N		N	
P87	PLASMA	N		N	
P88	STOOL	N		N	
P89	STOOL	N		N	
P90	STOOL	N		N	
P91	SPUTUM	N		N	
P92	SPUTUM	N		N	
Q1	SERUM	N		N	
Q2	STOOL	N		N	
Q3	STOOL	N		N	
Q4	STOOL	INC		N	
Q5	STOOL	N		N	
Q6	STOOL	N		N	
Q7	STOOL	N		N	
Q8	STOOL	N		N	
Q9	SERUM	N		N	
Q10	PLASMA	N		N	
Q11	STOOL	N		N	
Q12	STOOL	575	768	N	
Q13	STOOL	68		N	
Q14	STOOL	INC		N	
Q15	STOOL	N		N	



IMCB ID	DESC	1ST ARTUS	2ND ARTUS	PS 2 IMCB	IMCB Note
Q16	STOOL	N	N		
Q17	STOOL	N	N		
Q18	STOOL	206	1,423	POS	a->g (4656)
Q19	PLASMA	N	N		
Q20	SERUM	N	N		
Q21	PLASMA	N	N		
Q22	STOOL	N	N	POS	a->g (4656)
Q23	STOOL	456	325	POS	a->g (4656)
Q24	STOOL	453	477	N	
Q25	STOOL	236	N	N	
Q26	STOOL	1,083	719	N	
Q27	STOOL	786	469	N	
Q28	Stool: 2 Extract	850		N	
Q29	Stool: 2 Extract	N		POS	a->g (4656)
Q30	Stool: 2 Extract	N		N	
Q31	Stool: 2 Extract	INC		N	
Q32	Stool: 2 Extract	N		POS	a->g (4656)
Q33	Stool: 2 Extract	N		N	
Q34	Stool	1,326	519	N	
Q35	Stool	153	277	POS	a->g (4656)
Q36	Serum	N		N	
Q37	Stool: 2 Extract	1,253		N	
Q38	Stool: 2 Extract	N		N	
Q39	QS-1	?		N	
Q40	RNA EXTRACT+ci	?		N	LowerBand (not SAI)
Q41	2nd extract	606		POS	a->g (4656)
Q42	2nd extract	N		N	
Q43	stool	193	N	POS	a->g (4656)
Q44	stool	776	887	N	
Q45	Old PBS	250		N, N, N	repeat x2 negative
Q46	New PBS	N		N, N, N	repeat x2 negative
Q47	EBI + RNA Ex Cor	500		POS	Urbani, NEA
Q48	PLASMA	N		N	
Q49	STOOL	N		N	
Q50	PLASMA	N		N	
Q51	SUPERNATANT	N		N	
Q52	CELL PELLET	N		N	
Q53	STOOL	N		N	
Q54	STOOL	N		N	
Q55	STOOL	N		POS	a->g (4656)

IMCB ID	DESC	1ST ARTUS	2ND ARTUS	PS 2 IMCB	IMCB Note
Q56	SERUM	N		N	
Q57	SPUTUM (SUP)	N		N	
Q58	SPUTUM (C.P)	N		POS	a->g (4656)
Q59	STOOL	N		N	
Q60	PLASMA	N		N	
Q61	STOOL	N		N	
Q62	SPUTUM (SUP)	N		POS	a->g (4656)
Q63	SPUTUM (C.P)	N		N	
Q64	STOOL	N		N	
Q65	STOOL	N		POS	a->g (4656)
Q66	STOOL	N		N	
Q67	STOOL	N		N	
Q68	PLASMA	N		N	
Q69	STOOL	N		N	
Q70	SPUTUM	N		POS	(seq Q58)
Q71	SPUTUM	N		N	
Q72	STOOL	N		N	
Q73	SPUTUM	N		N	
Q74	SPUTUM	INC		N	
Q75	STOOL	N		N	
Q76	INHIBITOR TEST	(POS)		POSPOSPOS	Urbani, NEA
Q77	INHIBITOR TEST	(POS)		POSPOSPOS	Urbani, NEA

#### Example 8: SARS Real-Time PCR Diagnostic Kit (RT-PCR)

A kit of this example is typically prepared to contain 50 or 100 reactions. This Real-Time Reverse Transcriptase-Polymerase Chain Reaction (RT-PCR) kit is optimized to detect the presence of Severe Acute Respiratory syndrome Coronavirus Ribonucleic Acid (SARS coronavirus RNA) in a biological sample. This kit is optimized for use with the Applied Biosystems Real-Time PCR, ABI Prism 7500, but may be used with other suitable detection platforms as described elsewhere herein.

The portion of the SARS coronavirus genome amplified by the kit lies in the proteinase region of the SARS coronavirus RNA. The kit is sufficiently sensitive to detect a few molecules of RNA in each RT-PCR reaction.

*Components*

- This kit of this Example consists of the following 4 tubes:
- Tube 1 : Reaction Mix (e.g. ABI cat. No. 4309169)
- Tube 2 : Enzyme Mix (e.g. ABI cat. No. 4309169)
- Tube 3 : Probe Mix (3  $\mu$ M upper primer, 3  $\mu$ M lower primer,  
5                   2  $\mu$ M probe in 20 mM Tris, 1 mM EDTA pH 8.2)
- Tube 4 : Positive Control (RNA transcripts of the  
                                  gene targeted by the primers)

*Protocol*

## 1. RT-PCR

- 10           The following reaction mix is prepared in a 96-well  
Optical plate:

<u>Tube No.</u>	<u>Description</u>	<u>Vol/Rxn</u>
Tube 1	Reaction Mix	25.0 $\mu$ l
15 Tube 2	Enzyme Mix	1.25 $\mu$ l
Tube 3	Probe Mix	5.0 $\mu$ l
Tube 4	Distilled Water	13.75 $\mu$ l
-	RNA Sample	5.0 $\mu$ l
	Total Volume	50.0 $\mu$ l

- 20 Caution should be taken to avoid contamination.

## 2. - Thermal cycling conditions:

<u>Step</u>	<u>Temp (°C)</u>	<u>Duration</u>	<u>No. of Cycle(s)</u>
25 1	48	30 mins	1
2	95	10 mins	1
3	95	15 secs	50
4	60	60 secs	1

- 30 The IMCB-3 primer set and probe are tested for their ability  
to detect SARS coronavirus using the Stratagene real-time

PCR system Mx3000P. The system is used according to the manufacturer's instructions on samples from infected patients. The samples are diluted several fold to a total of viral copy number per 5 µl ranging between 7.5 to 6.

5 Results are shown in Table 4. The number of viral copies of SARS coronavirus per 5 µl ranges from 7.5 to 6 in all three runs. Row A (A1-3) is a control no virus sample. Row C (C1-C10) is detection of samples at 7.5 viral copies 5 µl, Row E (E1-E10) is detection of samples at 7.5 viral  
10 copies 5 µl, and Row G (G1-G10) is detection of samples at 7.5 viral copies 5 µl. The 4<sup>th</sup> column indicates how many positive results are detected per number of samples tested. The results demonstrate that the IMCB-3 primer set and probe provide a sensitive and specific assay for SARS coronavirus  
15 that is useful in a clinical setting.

Table 4: Sensitivity of the IMCB-3 primer set in real-time PCR analysis

Run #1:

Well	Well Name	Well Type	Threshold (dR)	Ct (dR)
A1	---	NTC	2907.809	0/3
A2	---	NTC	2907.809	
A3	---	NTC	2907.809	
C1	BMRC -7.5	Unknown	2907.809	8/10
C2	BMRC -7.5	Unknown	2907.809	
C3	BMRC -7.5	Unknown	2907.809	
C4	BMRC -7.5	Unknown	2907.809	
C5	BMRC -7.5	Unknown	2907.809	
C6	BMRC -7.5	Unknown	2907.809	
C7	BMRC -7.5	Unknown	2907.809	
C8	BMRC -7.5	Unknown	2907.809	
C9	BMRC -7.5	Unknown	2907.809	
C10	BMRC -7.5	Unknown	2907.809	
E1	BMRC -7.0	Unknown	2907.809	10/10
E2	BMRC -7.0	Unknown	2907.809	
E3	BMRC -7.0	Unknown	2907.809	
E4	BMRC -7.0	Unknown	2907.809	
E5	BMRC -7.0	Unknown	2907.809	
E6	BMRC -7.0	Unknown	2907.809	
E7	BMRC -7.0	Unknown	2907.809	
E8	BMRC -7.0	Unknown	2907.809	
E9	BMRC -7.0	Unknown	2907.809	
E10	BMRC -7.0	Unknown	2907.809	
G1	BMRC -6.5	Unknown	2907.809	10/10
G2	BMRC -6.5	Unknown	2907.809	
G3	BMRC -6.5	Unknown	2907.809	
G4	BMRC -6.5	Unknown	2907.809	
G5	BMRC -6.5	Unknown	2907.809	
G6	BMRC -6.5	Unknown	2907.809	
G7	BMRC -6.5	Unknown	2907.809	
G8	BMRC -6.5	Unknown	2907.809	
G9	BMRC -6.5	Unknown	2907.809	
G10	BMRC -6.5	Unknown	2907.809	

Run #2

report 040902

Well	Well Name	Well Type	Threshold (dr)	Ct (dr)
A1	NTC	2929.155	No Ct	
A2	NTC	2929.155	No Ct	
A3	NTC	2929.155	No Ct	0/3
C1	BMRC -7.5	Unknown 2929.155	41.05	
C2	BMRC -7.5	Unknown 2929.155	41	
C3	BMRC -7.5	Unknown 2929.155	42.06	
C4	BMRC -7.5	Unknown 2929.155	42.32	
C5	BMRC -7.5	Unknown 2929.155	42.22	9/10
C6	BMRC -7.5	Unknown 2929.155	41.92	
C7	BMRC -7.5	Unknown 2929.155	41.45	
C8	BMRC -7.5	Unknown 2929.155	40.72	
C9	BMRC -7.5	Unknown 2929.155	No Ct	
C10	BMRC -7.5	Unknown 2929.155	42.39	
E1	BMRC -7.0	Unknown 2929.155	42.03	
E2	BMRC -7.0	Unknown 2929.155	41.15	
E3	BMRC -7.0	Unknown 2929.155	41.11	
E4	BMRC -7.0	Unknown 2929.155	39.2	
E5	BMRC -7.0	Unknown 2929.155	No Ct	9/10
E6	BMRC -7.0	Unknown 2929.155	39.58	
E7	BMRC -7.0	Unknown 2929.155	43.15	
E8	BMRC -7.0	Unknown 2929.155	39.61	
E9	BMRC -7.0	Unknown 2929.155	40.38	
E10	BMRC -7.0	Unknown 2929.155	40.86	
G1	BMRC -6.5	Unknown 2929.155	38.57	
G2	BMRC -6.5	Unknown 2929.155	39.11	
G3	BMRC -6.5	Unknown 2929.155	38.63	
G4	BMRC -6.5	Unknown 2929.155	38.14	
G5	BMRC -6.5	Unknown 2929.155	37.77	10/10
G6	BMRC -6.5	Unknown 2929.155	38.33	
G7	BMRC -6.5	Unknown 2929.155	39.1	
G8	BMRC -6.5	Unknown 2929.155	38.08	
G9	BMRC -6.5	Unknown 2929.155	39.07	
G10	BMRC -6.5	Unknown 2929.155	39.04	

Run #3

report 040903

Well	Well Name	Well Type	Threshold (dr)	Ct (dr)
A1	NTC	3962.797	No Ct	
A2	NTC	3962.797	No Ct	
A3	NTC	3962.797	No Ct	0/3
C1	BMRC -7.5	Unknown 3962.797	43.21	
C2	BMRC -7.5	Unknown 3962.797	40.46	
C3	BMRC -7.5	Unknown 3962.797	42.1	
C4	BMRC -7.5	Unknown 3962.797	No Ct	
C5	BMRC -7.5	Unknown 3962.797	41.78	
C6	BMRC -7.5	Unknown 3962.797	40.7	5/10
C7	BMRC -7.5	Unknown 3962.797	No Ct	
C8	BMRC -7.5	Unknown 3962.797	No Ct	
C9	BMRC -7.5	Unknown 3962.797	No Ct	
C10	BMRC -7.5	Unknown 3962.797	No Ct	
E1	BMRC -7.0	Unknown 3962.797	38.4	
E2	BMRC -7.0	Unknown 3962.797	40.14	
E3	BMRC -7.0	Unknown 3962.797	40.03	
E4	BMRC -7.0	Unknown 3962.797	39.04	
E5	BMRC -7.0	Unknown 3962.797	41.12	10/10
E6	BMRC -7.0	Unknown 3962.797	39.48	
E7	BMRC -7.0	Unknown 3962.797	39.55	
E8	BMRC -7.0	Unknown 3962.797	39.89	
E9	BMRC -7.0	Unknown 3962.797	39.57	
E10	BMRC -7.0	Unknown 3962.797	42.92	
G1	BMRC -6.5	Unknown 3962.797	38.41	
G2	BMRC -6.5	Unknown 3962.797	37.41	
G3	BMRC -6.5	Unknown 3962.797	38.12	
G4	BMRC -6.5	Unknown 3962.797	37.73	
G5	BMRC -6.5	Unknown 3962.797	37.6	10/10
G6	BMRC -6.5	Unknown 3962.797	38.58	
G7	BMRC -6.5	Unknown 3962.797	37.9	
G8	BMRC -6.5	Unknown 3962.797	38.38	
G9	BMRC -6.5	Unknown 3962.797	38.27	
G10	BMRC -6.5	Unknown 3962.797	38.1	

It will be understood by those of skill in the art that the presence or absence of components in the above method, the concentrations of chemicals, the cycling conditions, and the equipment can be modified to suit the particular needs and to optimize reaction conditions.

Moreover, those of skill in the art will recognize that the above method and the above description can encompass modifications that fall within the spirit and scope of the instant invention.

#### References

The following references are cited in the instant specification. Each of the following references is hereby incorporated by reference in its entirety and for all purposes by such citation:

Drosten, C., et al., 2003. Identification of a novel coronavirus in patients with severe acute respiratory syndrome. N. Engl. J. Med. 348:1967-1976.

Notomi, T. et al., 2000. Loop-mediated isothermal amplification of DNA. Nucleic Acids Research 15:E63.

Rota, P. A., et al. 2003. Characterization of a novel coronavirus associated with severe acute respiratory syndrome. Science 300:1394-1399.

Ruan, Y.J., et al., 2003. Comparative full-length genome sequence analysis of 14 SARS coronavirus isolates and common mutations associated with putative origins of infection. Lancet. 361, 1779-1785.

WHO Update 71. Status of diagnostic tests, training course  
in China. [http://www.who.int/entity/csr/don/2003\\_06\\_02\\_a/en](http://www.who.int/entity/csr/don/2003_06_02_a/en)

**This Page is Inserted by IFW Indexing and Scanning  
Operations and is not part of the Official Record**

**BEST AVAILABLE IMAGES**

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images include but are not limited to the items checked:

- ☐ BLACK BORDERS
- ☐ IMAGE CUT OFF AT TOP, BOTTOM OR SIDES
- ☐ FADED TEXT OR DRAWING
- ☒ BLURRED OR ILLEGIBLE TEXT OR DRAWING
- ☐ SKEWED/SLANTED IMAGES
- ☐ COLOR OR BLACK AND WHITE PHOTOGRAPHS
- ☐ GRAY SCALE DOCUMENTS
- ☐ LINES OR MARKS ON ORIGINAL DOCUMENT
- ☐ REFERENCE(S) OR EXHIBIT(S) SUBMITTED ARE POOR QUALITY
- ☐ OTHER: \_\_\_\_\_

**IMAGES ARE BEST AVAILABLE COPY.**

**As rescanning these documents will not correct the image problems checked, please do not report these problems to the IFW Image Problem Mailbox.**